



Blood Plasma

1. Freeze (at -20C) plasma samples immediately after separation from other blood components (plasma top layer; leucocytes & erythrocytes bottom layers)
2. Transfer 500 μ L of sample to a 7-mL test tube and freeze-dry over night
3. Add 5 mL of petroleum ether to each test tube and cap tightly
4. Sonicate with heat "on" samples under a fume hood for 15 minutes
5. Decant the used solution
6. Repeat steps 3-5 for a second and third rinse
7. Add 5mL of deionized to each test tube and cap tightly
8. Sonicate with heat "on" samples under a fume hood for 15 minutes
9. Decant the used water
10. Repeat steps 7-9 two more times for a second and third water rinse
11. Freeze-dry samples over night to remove any remaining water
12. Homogenize the extracted tissue within the test tube and weigh ~1.0 mg of the sample into packing tin
13. Combust at 1100C and deliver gases via continuous-flow for analysis of $\delta^{13}C$, $\delta^{15}N$, percentage carbon and percentage nitrogen using and isotope ratio mass spectrometer (IRMS).
14. Calculate atomic C:N ratio from using the measured mass of each element:
 - a. $C:N = (\text{mass of C})/(\text{mass of N})$
 - b. If C:N is greater than 3.0, it likely still contains lipids and another aliquot of sample should be extracted using additional petroleum ether rinses
 - c. Report mean C:N and standard errors for tissues subject to SIA



Whole Blood

1. Freeze (at -20C) samples immediately after separation from other blood components (plasma top layer; leucocytes & erythrocytes bottom layers)
2. Transfer 250 μ L of sample to a 7-mL test tube and freeze-dry over night
3. Add 5 mL of petroleum ether to each test tube and cap tightly
4. Sonicate with heat "on" samples under a fume hood for 15 minutes
5. Decant the used solution
6. Repeat steps 3-5 for a second and third rinse
7. Add 5mL of deionized to each test tube and cap tightly
8. Sonicate with heat "on" samples under a fume hood for 15 minutes
9. Decant the used water
10. Repeat steps 7-9 two more times for a second and third water rinse
11. Freeze-dry samples over night to remove any remaining water
12. Homogenize the extracted tissue within the test tube and weigh ~1.0 mg of the sample into packing tin
13. Combust at 1100C and deliver gases via continuous-flow for analysis of $\delta^{13}\text{C}$, $\delta^{15}\text{N}$, percentage carbon and percentage nitrogen using and isotope ratio mass spectrometer (IRMS).
14. Calculate atomic C:N ratio from using the measured mass of each element:
 - a. $\text{C:N} = (\text{mass of C})/(\text{mass of N})$
 - b. If C:N is greater than 3.0, it likely still contains lipids and another aliquot of sample should be extracted using additional petroleum ether rinses
 - c. Report mean C:N and standard errors for tissues subject to SIA



Muscle/Skin Tissue

1. Freeze (at -20C) samples immediately after collection
2. Freeze-dry whole sample in original containers for ~48 hrs
3. Grind whole sample with ball mill for 3 min
4. Transfer >5 mg (dry weight) of ground sample into 15 mL test tubes with PE caps
5. Add 5 mL of petroleum ether to each test tube and cap tightly
6. Sonicate with heat "on" samples under a fume hood for 15 minutes
7. Decant the used solution
8. Repeat steps 3-5 for a second and third rinse
9. Add 5mL of deionized to each test tube and cap tightly
10. Sonicate with heat "on" samples under a fume hood for 15 minutes
11. Centrifuge all tubes for 2min
12. Decant the used water
13. Repeat steps 7-9 two more times for a second and third water rinse
14. Freeze-dry samples over night to remove any remaining water
15. Homogenize the extracted tissue within the test tube and weigh ~1.0 mg of the sample into packing tin
16. Combust at 1100C and deliver gases via continuous-flow for analysis of $\delta^{13}C$, $\delta^{15}N$, percentage carbon and percentage nitrogen using and isotope ratio mass spectrometer (IRMS).
17. Calculate atomic C:N ratio from using the measured mass of each element:
 - a. $C:N = (\text{mass of C})/(\text{mass of N})$
 - b. If C:N is greater than 3.0, it likely still contains lipids and another aliquot of sample should be extracted using additional petroleum ether rinses
 - c. Report mean C:N and standard errors for tissues subject to SIA

References:

Hussey et al 2012 Stable isotopes and elasmobranchs: tissue types, methods, applications and assumptions. *Journal of Fish Biology* (2012) 80, 1449–1484
doi:10.1111/j.1095-8649.2012.03251.x

Kim, S. L. & Koch, P. L. (2012). Methods to collect, preserve, and prepare elasmobranch tissues for stable isotope analysis. *Environmental Biology of Fishes* 95:53-63. doi: 10.1007/s10641-011-9860-9

Kim et al 2012 Carbon and nitrogen discrimination factors for elasmobranch soft tissues based on a long-term controlled feeding study. *Environ Biol Fish* (2012) 95:37–52 DOI 10.1007/s10641-011-9919-7

Newsome SD, Clementz MT, Koch PL (2010) Using stable isotope biogeochemistry to study marine mammal ecology. *Mar Mamm Sci* 26:509–572