



Protocol for Soil Extraction of Nitrate-N and Ammonium-N

Equipment

- 1) 100 ml Erlenmeyer flasks & clean rubber stoppers
- 2) Whatman #42 filter paper
- 3) Long-stemmed funnels and funnel rack
- 4) Mechanical shaker box
- 5) Polyethylene scintillation vials with poly-lined caps

2M potassium chloride extraction solution

- 1) in 1000 ml volumetric flask, add
150 g KCl
800 ml deionized water (diH₂O)
- 2) Add stir bar; place on stir plate and mix until dissolved.
- 3) Remove stir bar and bring to full 1000 ml volume with diH₂O.
- 4) Cap with parafilm and invert several times to mix.

Procedure

- 1) Place 4 g wet-weight soil into each flask; add 20 ml extraction solution to each. Stopper flasks.
- 2) Place flasks in shaker box and shake 1 hour, medium speed.
- 3) Filter flask contents through funnels lined with filter paper; catch filtrate in labelled scintillation vials.
- 4) Cap vials and refrigerate until ready to analyze.
- 5) Analyze extracts by (Alpkem) continuous-flow colorimetry.

Post-analysis calculations

Calculate soil content from the extract's nitrate and ammonium values as determined by colorimetric analysis:

$$\text{ug analyte per g dry soil} = (\text{determined value, mg/L}) \times \frac{1000 \text{ ug}}{1 \text{ mg}} \times \frac{0.02 \text{ L}}{\text{soil dry weight (gm)}}$$

(The dry weight value of the soil sample is based on percent moisture data for each sample, obtained separately.)

Bibliography

Keeney, D. R. and D. W. Nelson. 1987.

Nitrogen--Inorganic Forms, sec. 33-3, extraction of exchangeable ammonium, nitrate, and nitrite. pp.648-9. In A. L. Page et al., eds., Methods of Soil Analysis: Part 2, Chemical and Microbiological Properties. Agronomy, A Series of Monographs, no.9 pt.2, Soil Science Society of America, Madison, Wisconsin USA.